

UNIMED
SUTURES

UNICRYL M

UNICRYL M

UNICRYL M

- Monofilament Polydioxanone E.P. Violet and Undyed.
- Synthetic absorbable sterile surgical suture prepared from polyester p-dioxanone.
- Available in range of gauge sizes and lengths, non-needled or attached to stainless steel needles of varying types and sizes.
- Minimal tissue reaction.
- Easy passage through tissue.
- Intended for use where an absorbable suture or ligature is indicated. It can be used in all types of tissue approximation, pediatric cardiovascular tissue and ophthalmic surgery.
- Sterilised by Ethylene Oxide.
- Sutures are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), European Pharmacopoeia (Ph. Eur.), and the United States Pharmacopoeia (U.S.P.).

UNICRYL M - Approximate Breaking Strength Retention

