

UNIMED
SUTURES

UNICHROM

UNICHROM

UNICHROM

- Chromic Catgut B.P.
- Manufactured from uniform, twisted strands of collagen which are treated to retain tensile strength longer than plain catgut.
- Packaged in a hydrating fluid containing isopropanol and water.
- Origin - purified animal intestines (bovine).
- Type - absorbable by enzymatic digestion.
- Sterilised by Gamma Irradiation.
- Sutures are guaranteed to meet the general requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur) and the United States Pharmacopoeia (U.S.P.).