









































































UNICHROM CHROMICISED CATGUT B.P.

NEEDLE			LENGTH	GAUGE						
				6/0 (1)	5/0 (1.5)	4/0 (2)	3/0 (3)	2/0 (3.5)	0 (4)	
13mm										
Curved round bodied			60cm			DW 381				
16mm										
1/2 circle round bodied			75cm	DW 815	DW 816	DW 817	DW 818	DW 819		
Curved cutting			75cm			DW 811	DW 812	DW 813		
Curved cutting double			38cm				DW 812D			
20mm										
1/2 circle round bodied			75cm	DW 559	DW 560	DW 561	DW 562	DW 553		
1/2 circle taper cut			75cm				DT 222	DT 223		
Curved cutting - fine			45cm			DW 658				
Curved reverse cutting - fine			45cm			DW 658R				
25mm										
1/2 circle round bodied			75cm			DW 501	DW 502	DW 503	DW 504	
Curved cutting - fine			45cm			DW 571				
1/2 circle taper cut			75cm				DT 592	DT 593	DT 594	
1/2 circle taper cut - heavy			75cm				DT 592H	DT 593H		















UNICHROM CHROMICISED CATGUT B.P.

NEEDLE			LENGTH	GAUGE					
				4/0 (2)	3/0 (3)	2/0 (3.5)	0 (4)	1 (5)	2 (6)
30mm									
Curved round bodied			75cm			DW 403	DW 404		
1/2 circle round bodied			75cm	DW 511	DW 512	DW 513	DW 514	DW 515	
1/2 circle round bodied - heavy			75cm				DW 1014	DW 1015	DW 1016
1/2 circle taper cut			75cm		DT 412		DT 414	DT 415	
1/2 circle taper cut - heavy			75cm					DW 415H	
1/2 circle cutting			75cm			DW 573			
1/2 circle trocar point - heavy			75cm					DW 1025	
35mm									
1/2 circle round bodied			75cm			DW 843	DW 844	DW 845	DW 846
1/2 circle round bodied - double			75cm			DW 843D			
1/2 circle round bodied - heavy			75cm				DW 844H	DW 845H	DW 846H
1/2 circle taper cut			75cm			DT 873	DT 874		
1/2 circle taper cut - heavy			75cm				DT 874H	DT 875H	DT 876H
1/2 circle trocar point - heavy			75cm					DW 895	
38mm									
1/2 circle round bodied			75cm			DW 883	DW 884	DW 885	
1/2 circle round bodied - heavy			75cm						DW 886H









UNICHROM CHROMICISED CATGUT B.P.

NEEDLE		LENGTH	GAUGE						
			4/0 (2)	3/0 (3)	2/0 (3.5)	0 (4)	1 (5)	2 (6)	
40mm									
1/2 circle round bodied			75cm		DW 532	DW 533	DW 534	DW 535	
1/2 circle round bodied - fine			75cm					DW 525F	
1/2 circle round bodied - heavy			75cm			DW 1033	DW 1034	DW 1035	DW 1036
1/2 circle reverse cutting			75cm				DW 734R		DW 736R
1/2 circle cutting			75cm		DW 732				
1/2 circle taper cut			75cm					DT 497	
1/2 circle taper cut - heavy			75cm					DT 1045	
1/2 circle taper cut - heavy double			60cm					DT 1045D	
1/2 circle Trocar point - heavy			75cm					DW 525	DW 526





UNICHROM CHROMICISED CATGUT B.P.

NEEDLE		LENGTH	GAUGE						
			4/0 (2)	3/0 (3)	2/0 (3.5)	0 (4)	1 (5)	2 (6)	
45mm									
Curved round bodied			75cm					DW 435	
1/2 circle round bodied			75cm			DW 943	DW 944		
1/2 circle round bodied - heavy			75cm			DW 953	DW 954	DW 955	DW 956
1/2 circle round bodied - heavy/multi			75cm					DW 955M	DW 956M
1/2 circle cutting			75cm				DW 974		
1/2 circle taper cut - heavy			75cm			DT 993	DT 994	DT 995	DT 996
1/2 circle Trocar point - heavy			75cm					DW 1065	DW 1066

UNICHROM CHROMICISED CATGUT B.P.

NEEDLE	LENGTH	GAUGE						
		4/0 (2)	3/0 (3)	2/0 (3.5)	0 (4)	1 (5)	2 (6)	
50mm								
1/2 circle round bodied 	 75cm				DW 563		DW 565	
1/2 circle round bodied - double 	 75cm				DW 563D		DW 566D	
1/2 circle round bodied heavy 	 75cm					DW 1054	DW 1055	DW 1056
1/2 circle cutting 	 75cm					DW 1064		

UNICHROM CHROMICISED CATGUT B.P.

NEEDLE	LENGTH	GAUGE	
		1 (5)	2 (6)
60mm			
1/2 circle round bodied 	 75cm		DW 586
63mm			
Curved round bodied - blunt 	 75cm	DW 1225	
65mm			
1/2 circle round bodied 	 75cm		DW 1266
85mm			
Curved round bodied - blunt 	 75cm	DW 1355	

NON-NEEDED UNICHROM	LENGTH	GAUGE							
		6/0 (1)	5/0 (1.5)	4/0 (2)	3/0 (3)	2/0 (3.5)	0 (4)	1 (5)	2 (6)
Ligature sets	5 x 38cm				DW 52	DW 53	DW 54	DW 55	DW 56
Cut lengths	75cm	DW 69	DW 70	DW 71	DW 72	DW 73	DW 74	DW 75	DW 76
Standard lengths	1.5m		DW 200	DW 201	DW 202	DW 203	DW 204	DW 205	DW 206