

**UNIMED**  
SUTURES

UNIESTER

**UNIESTER**

# UNIESTER

---

- Braided polyster suture B.P., White.
- Braided polyster suture B.P., Black.
  
- Type - non-absorbable synthetic material.
  
- Impregnated to reduce capillary action and to minimise tissue trauma.
  
- High tensile strength for reliability and support.
  
- Unaffected by wetting or implantation in tissue.
  
- The flexibility of a braided construction.
  
- Uses - many cardiovascular, vascular and orthopaedic procedures where a strong, braided suture is indicated.
  
- Sterilised by gamma Irradiation.
  
- Sutures are guaranteed to meet the general requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur) and the United States Pharmacopoeia (U.S.P.).