

UNIMED
SUTURES

UNIESTER C

UNIESTER C

UNIESTER C

- Coated Braided Polyester Suture B.P., Green
- Manufactured from a smooth, braided polyester material with a unique coating.
- Type - non-absorbable synthetic material.
- Extremely strong for long term support.
- Soft and pliable providing good knot security.
- Minimal tissue reaction.
- Easy passage through tissue and reduced tissue trauma.
- Uses - cardiac,, vascular and orthopaedic procedures where long term support from a braided non - absorbable suture is essential.
- Sterilised by gamma Irradiation.
- C sutures are guaranteed to meet the general requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur) and the United States Pharmacopoeia (U.S.P.).