

**UNIMED**  
SUTURES

UNIPLAIN

***UNIPLAIN***

# UNIPLAIN

---

- Catgut Plain B.P.
- Manufactured from uniform, twisted strands of collagen.
- Packaged in a hydrating fluid containing isopropanol and water.
- Origin - purified animal intestines (bovine).
- Type - absorbable by enzymatic digestion.
- Uses - UNIPLAIN sutures are intended for use where an absorbable suture or ligature is indicated.
- Sterilised by Gamma Irradiation.
- Sutures are guaranteed to meet the general requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur) and the United States Pharmacopoeia (U.S.P.).