

UNIMED
SUTURES

UNIPRO

UNIPRO

UNIPRO

- Monofilament Polypropylene suture B.P. Blue.
- Synthetic non-absorbable sterile surgical suture.
- Available in range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes.
- Inert - Minimal tissue reaction.
- Very smooth for easy passage through tissue.
- Less memory, very pliable-good handling.
- Easy knot tying, excellent knot security.
- Even in infected tissue, polypropylene provides long term strength retention.
- Intended for use in procedures where a smooth monofilament suture offering very long term support is required. It can be used in general soft tissue approximation or ligation. Also, recommended for cardiovascular, ophthalmic and neurological procedures.
- Sterilised by Ethylene Oxide.
- Sutures are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur) and the United States Pharmacopoeia (U.S.P.).