

UNIMED
SUTURES

UNISTEEL

UNISTEEL

UNISTEEL

- Multifilament Braided Stainless Steel Wire Suture Ph. Eur.
- Monofilament Stainless Steel Wire Suture Ph. Eur.

- Type - non-absorbable.

- Both braided and monofilament stainless steel wire sutures comply to BS7252.

- Monofilament - very strong for sternal application.
powerful cutting needle for smooth passage through bone.

- Multifilament - insulated, soft braided wire for flexibility.

- Monofilament - sternal closure.
multifilament - temporary pacemaker wires.

- Sterilised by Gamma Irradiation.

- Sutures are guaranteed to meet the general requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur) and the United States Pharmacopoeia (U.S.P.).