

UNIMED
SUTURES

UNIVYL

UNIVYL

UNIVYL

- Monofilament Polyvinylidene Fluoride Suture.
- Type - non-absorbable synthetic material.
- Excellent typing properties compared with other monofilaments.
- Very strong
- Less memory than other monofilaments.
- Intended for use wherever a strong, smooth monofilament suture, offering long term support is indicated.
- Sterilised by Gamma Irradiation.
- Sutures are guaranteed to meet the general requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph. Eur) and the United States Pharmacopoeia (U.S.P.).